

MDS-Vet™ offers two options for your practice: the Exam Room Otoscope System and the Video Otoscope System featuring Auri-Kleen™. Both otoscope systems provide convenient, easy-to-use and reliable solutions for veterinary practices from the MDS-Vet™ line of products. These systems provide the end-user with a versatile and valuable tool to perform both Otoscopy and Rhinoscopy with a single endoscope.


EXAM ROOM OTOSCOPE SYSTEM - ITEM # 48008



FEATURES:

- High quality optical endoscope
- Compact otoscopy system with Color CCD Camera
- Compatible with a laptop or tablet via USB port
- Converts to 2.7 mm rhinoscope
- 50-Watt LED Light Source and cable (110V-240V)
- Portable LED Light Source included with every scope purchase


VIDEO OTOSCOPE SYSTEM WITH AURI-KLEEN™ - ITEM # 48055




FEATURES:

- Includes everything in Exam Room Otoscope System
- Adds the Auri-Kleen™ Irrigation / Aspiration Machine
- Pair with optional cart for easy storage and transport
- Pretreat with OtoKleen™ Ear Cleansing Solution to help lift dirt and debris before deep ear flush
- Pair with OtoKleen™ Ear Cleaning Cannulas for direct delivery into the ear canal - luer locks onto the Auri-Kleen™ system

Item	Product Description	UOM
48008	Exam Room Video Otoscope System	Each
48055	Video Otoscope System with Auri-Kleen™ Irrigation / Aspiration Machine	Each
48100	Auri-Kleen™ Irrigation / Aspiration Machine	Each
48009	2.7 mm x 105 mm Rhinoscope	Each
48010	Otoscope Conversion Sheath	Each
48405	Alligator Foreign Body Forceps, 4.5FR, 8" Length	Each
48304	50-Watt LED Light Source	Each
48217	Light Guide for 50-Watt Light Source	Each
48050	Cart, Welded Metal with Drawer (optional add-on)	Each
48201	High-Resolution USB Endoscopy Camera	Each
48203	Focusing Endoscope Coupler	Each
48335	Portable LED Light Source	Each
48104	OtoKleen™ Ear Cleaning Cannulas, Small 8FR	10/Bag
48103	OtoKleen™ Ear Cleansing Solution - 16 oz	Each



Before



After

**EVERYTHING YOU
NEED TO CREATE A
PROFIT CENTER IN
YOUR PRACTICE.**



Deep ear flushing made easy utilizing an integrated irrigation and suction device with hands-free operation.

UNIQUE, HANDS-FREE, FOOT OPERATED CONTROL


- Dual function foot pedal provides irrigation and aspiration for simplified deep ear flushing
- High capacity (1.5 Liter) bottle with optimum irrigation pressure for dispensing fluid into patient's ear
- Medical grade, durable suction reservoir collects large amounts of aural debris and waste; easy to clean and sterilizable

SPECIFICALLY DESIGNED FOR VETERINARY / PROFESSIONAL USAGE

- Designed specifically for a single operator
- Ideal system to complement a handheld or video otoscope
- Air cooled system for extended operator usage
- Dimensions: 13" x 16" x 13" (with bottle)



Extreme wax build-up before flushing



Cleaning ear canal and tympanic membrane using Auri-Kleen™

Item	Product Description	UOM
48100	Auri-Kleen™ Irrigation / Aspiration Machine	Each
48105	Replacement Suction Canister, 1,500 ml	Each
48107	Replacement Silicone Tubing Assembly for Auri-Kleen™ Irrigation / Aspiration Machine, 7 ft	Each

**FOOT PEDAL
OPERATED DESIGN!**